BIOLOY®

SURFACE

permedica
DENTAL

BIOLOY®

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
BIOLOY® is a titanium niobium nitride (TiNbN) coating applied by Physical Vapour Deposition (PVD).

TiNbN is a thin ceramic, monolayer, composite coating, extremely hard and with high resistance against wear, scratches and corrosion.

TiNbN acts as a surface coat to "hide" the titanium implant beneath.

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BIOLOY® coating is performed by Permedica S.p.A.

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
Thickness: $4,5 \pm 1,5 \mu\text{m}$ (Ball-Crater test - UNI EN 1071-2:2003).

Biocompatibility: Non cytotoxic (ISO 10993-5).

Roughness: $S_a = 0.7 \mu\text{m} \pm 0.04$. Surface roughness depends on substrate microtopography.

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X-Ray Photoelectron Spectroscopy spectrum of surface chemical composition of BIOLOY®.

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Permedica started to use BIOLOY® coating in orthopedic field in 2007 to coat joint prosthetic components.

The ten-year experience in clinical practice with Bioloy®-coated knee replacement prostheses by Permedica and the excellent clinical results obtained, give proven excellent performances and reliability of BIOLOY® coating.

TiNbN coating has equivalent characteristics and properties compared to titanium nitride (TiN), well-known PVD coating for biomedical applications. TiNbN coating differs from TiN by the chemical element niobium (Nb) which give higher hardness and density to the coating [1].


Some researches on surface-modified biomaterials recently addressed the osseointegration behavior of TiN coated-titanium surfaces, with in vitro and in vivo comparable [2-4] or better [5-7] results compared to uncoated titanium surfaces of similar roughness.

Moreover, TiN-coated titanium surfaces seem to prevent bacterial colonization by inhibition of microbial adhesion and proliferation. Some in vitro and in vivo studies showed how TiN-coated titanium surfaces have the ability to reduce oral bacteria adhesion and proliferation in comparison of uncoated titanium surface of comparable roughness [8-13].

Latest researches on the TiN-coated titanium surface ability to osseointegrate and its ability to reduce oral bacterial adhesion and biofilm formation, have suggested the use of TiNbN PVD coating with dental implants and, thus, the development of BIOLOY® Surface by Permedica trying to reduce periimplantitis risk.

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Titanium surface is micro-structured by grit-blasting with a calcium magnesium carbonate $[\text{CaMg}(\text{CO}_3)_2]$ based media and then PVD coated by a thin TiNbN layer.

Surface average roughness (S_a) of BIOLOY® Surface is $0.7 \mu\text{m}$, thus a moderately roughened surface micro topography.

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SEM image at higher magnification of the microtopography of the calcium magnesium carbonate blasted titanium surface before PVD coating.

Surface average roughness (S_a) is $0.8 \mu\text{m}$.

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
SEM image at higher magnification of the microtopography of the BIOLOY® Surface, after PVD coating: macro-particles, or droplets, fused together are visible on the coating surface.

The coating modifies the typical sandblasted surface topography by rounding peaks and valleys and by the presence of said droplets, but preserving at the same time the surface average roughness.

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

PVD is a coating method by which positively charged metal ions are evaporated in a vacuum chamber and react with nitrogen gas introduced to the chamber. The surfaces that are to be coated are negatively charged to allow a strong bond to form between the substrate (implants) and the coating.

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Evaporation of the coating metal from the solid state to the plasma state is carried out by a high power electric arc between the surface of the titanium-niobium cathode and the coating chamber (anode).

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Analysis by SEM imaging of adhesion and proliferation at 24, 48, 72 hours of osteoblastic-like SaOS-2 cells cultured onto titanium specimens treated with three different surface finishing methods: machined, sandblasted, sandblasted and then TiNbN coated (BIOLOY® Surface). Cells viability with Alamar Blue assay [14].

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


Machined titanium

Cells seeded onto machined titanium smooth surface, took a widely spread shape. Cytoplasmatic pseudopods are well visible. Cells converge and join together forming interconnecting bridges to stabilize on the surface (red arrows).

Sandblasted titanium

Cells seeded onto roughened surface of sandblasted titanium took a coiled up, enclosed shape over the surface peaks, spreading along the surface depressions (red arrows).

BIOLOY® Surface treated titanium


Cells seeded onto TiNbN-coated sandblasted titanium (BIOLOY® Surface) took a morphology characterized by both shapes seen on machined titanium and on sandblasted titanium surfaces. Cells spread over the smoother peaks and droplets of the BIOLOY® Surface.

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Evaluation of cell proliferation by SEM imaging at 24, 48, 72 hours

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Cell proliferation on specimen surface of SaOS-2 cell seeded and cultured for 24, 48, 72 hours (magnification 550x).

Red areas are surface portions of the specimen uncovered by cells.

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
Cells proliferate progressively covering the specimen surfaces. At 72 hours the surfaces were completely covered by cells in both machined titanium and BIOLOY® Surface.

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Cell viability has been tested by AlamarBlue® assay. AlamarBlue® is a non-cytotoxic redox indicator that produces a color change from blue into red (fluorescence signal) in response to cell metabolic activity and proliferation.

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
*Cell viability of cells cultured on TiNbN-coated sandblasted titanium (BIOLOY® Surface) is statistically higher than cell viability of cells cultured on both sandblasted titanium and machined titanium at 48 and 72 hours (*p-value < 0.05; **p-value < 0.001).*

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
TiNbN coating (BIOLOY® Surface) promotes cell confluence and bonding by reducing peaks from sandblasted surface and, at the same time, provides fixation spots necessary for cell adhesion. This surface characterization increases osteoblastic-like cell viability and proliferation.

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Magnified image of a sample stained with Toluidine Blue and Pyronin Y which shows numerous newly formed blood vessels inside a new formed osteon (red arrows). Total magnification 400x.

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Magnified image of a sample stained with Toluidine Blue and Pyronin Y which shows in detail a soft tissue part of the sample. Note the physiological presence of some inflammatory cells such as granulocytes (red arrows), but the presence of an inflammatory infiltrate can be excluded. Total magnification 400x.

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The histological aspect of the regenerated periimplant bone around TiNbN-coated (BIOLOY® Surface) mini-implants has showed an active remodeling with a direct contact to the implant surface:

- *Bone is highly organized in lamellae with high degree of mineralization.*
- *Medullary spaces are rich in new blood vessels without inflammatory infiltrate.*

TiNbN-coated (BIOLOY® Surface) dental implants have showed favorable results in terms of biocompatibility and osseointegration.

BIOLOY® SURFACE-TREATED DENTAL IMPLANTS

«BIOLOY® Surface has the ability to osseointegrate as well as roughened titanium surfaces»

«Potential ability of BIOLOY® Surface to reduce oral bacteria adhesion and biofilm formation»

permedica
DENTAL

BIBLIOGRAFIA

REFERENCES

1. Serro AP, Completo C, Colaço R, dos Santos F, Lobato da Silva C, Cabral JMS, Araújo H, Pires E, Saramago B. A comparative study of titanium nitrides, TiN, TiNbN and TiCN, as coatings for biomedical applications. *Surface and Coatings Technology* 2009 Sep;203,(24):3701–07.
2. Scarano A, Piattelli M, Vrespa G, Petrone G, Iezzi G, Piattelli A. Bone healing around titanium and titanium nitride-coated dental implants with three surfaces: an experimental study in rats. *Clin Implant Dent Relat Res*. 2003;5(2):103-11.
3. Ferraz EP, Sa JC, de Oliveira PT, Alves C Jr, Beloti MM, Rosa AL. The effect of plasma-nitrided titanium surfaces on osteoblastic cell adhesion, proliferation, and differentiation. *J Biomed Mater Res A*. 2014 Apr;102(4):991-8.
4. Ferraz EP, Sverzut AT, Freitas GP, Sá JC, Alves C Jr, Beloti MM, Rosa AL. Bone tissue response to plasma-nitrided titanium implant surfaces. *J Appl Oral Sci*. 2015 Jan-Feb;23(1):9-13.
5. Durual S, Pernet F, Rieder P, Mekki M, Cattani-Lorente M, Wiskott HW. Titanium nitride oxide coating on rough titanium stimulates the proliferation of human primary osteoblasts. *Clin Oral Implants Res*. 2011 May;22(5):552-9.
6. Durual S, Rieder P, Garavaglia G, Filieri A, Cattani-Lorente M, Scherrer SS, Wiskott HW. TiNOx coatings on roughened titanium and CoCr alloy accelerate early osseointegration of dental implants in minipigs. *Bone*. 2013 Jan;52(1):230-7.
7. Sovak G, Gotman I, Weiss A. Osseointegration of Ti-6Al-4V alloy implants with a titanium nitride coating produced by a PIRAC nitriding technique: a long-term time course study in the rat. *Microsc Microanal*. 2015 Feb;21(1):179-89.
8. Ji MK, Park SW, Lee K, Kang IC, Yun KD, Kim HS, Lim HP. Evaluation of antibacterial activity and osteoblast-like cell viability of TiN, ZrN and (Ti1-xZrx)N coating on titanium. *J Adv Prosthodont*. 2015 Apr;7(2):166-71.
9. Velioglu N, Akova T, Ozkomur A. Effects of Hard Thin-Film Coatings on Adhesion of Early Colonizer Bacteria Over Titanium Surfaces. *Implant Dent*. 2016 Feb;25(1):114-21.
10. Annunziata M, Oliva A, Basile MA, Giordano M, Mazzola N, Rizzo A, Lanza A, Guida L. The effects of titanium nitride-coating on the topographic and biological features of TPS implant surfaces. *J Dent*. 2011 Nov;39(11):720-8.
11. Groessner-Schreiber B, Hannig M, Dück A, Griepentrog M, Wenderoth DF. Do different implant surfaces exposed in the oral cavity of humans show different biofilm compositions and activities? *Eur J Oral Sci*. 2004 Dec;112(6):516-22.
12. Grössner-Schreiber B, Griepentrog M, Hausteil I, Müller WD, Lange KP, Briedigkeit H, Göbel UB. Plaque formation on surface modified dental implants. An in vitro study. *Clin Oral Implants Res*. 2001 Dec;12(6):543-51.
13. Scarano A, Piattelli M, Vrespa G, Caputi S, Piattelli A. Bacterial adhesion on titanium nitride-coated and uncoated implants: an in vivo human study. *J Oral Implantol*. 2003;29(2):80-5.
14. Study report. Department of Biomedical, Surgical and Dental Sciences, Università degli Studi di Milano, Milan, Italy.
15. Frattini D, Canciani E, Dellavia C. Biocompatibility and histologic evaluation of the osseointegration of dental implants coating with titanium niobium nitride. *Clin. Oral Impl. Res*. 2016;27 (Suppl. 13).

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permedica deutschland GmbH