EXACTIVE®

DENTAL IMPLANT SYSTEM

permedica
DENTAL

EXACTIVE®

INDICE

INDEX

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

BIOLOY®

ELEVATA ESTETICA
HIGH AESTHETIC

CONNESSIONE CONICA
CONICAL CONNECTION

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GUIDELINES FOR SELECTING THE IMPLANT

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
In implantology treatments, it is preferable to use implants whose diameter is proportionate to the size of the missing element, optimising the aesthetic and biomechanical result.

EXACTIVE® implants can be used for any quadrant, especially in bone classes with density D1-D2 and also in bone classes D3-D4 for their stabilisation ability that is higher in cases of reduced primary stability due to bone density.

In the following table, according to the tooth to be replaced, the diameters of the implants are recommended and shown:

Mascella Superiore Upper Jaw	Ø Impianti EXACTIVE® EXACTIVE® Implant
Incisivi Centrali / Central Incisors	3.85 / 4.3
Incisivi Laterali / Lateral Incisors	3.3 / 3.85 / 4.3
Canini / Canines	3.85 / 4.3 / 5.0
Premolari / Premolars	3.85 / 4.3 / 5.0
Molari / Molars	4.3 / 5.0 / 6.0
Mandibola Raw	Ø Impianti EXACTIVE® EXACTIVE® Implant
Incisivi Centrali / Central incisors	3.3 / 3.85
Incisivi Laterali / Lateral incisors	3.3 / 3.85
Canini / Canines	3.85 / 4.3 / 5.0
Premolari / Premolars	3.85 / 4.3 / 5.0
Molari / Molars	4.3 / 5.0 / 6.0

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

permedica is one of the few companies in Italy that uses decontamination procedures of a high technological and qualitative level, which, in addition to a high degree of cleanliness, also guarantee the high wettability of the fixture documented both by the immediate "inks tests" as well as by specific tests to measure the angle of contact, which are carried out at our internal laboratory.

The plasma decontamination process of argon allows for an effective decontamination and generates the maximum degree of surface cleaning.

The plasma, generated inside a reactor, consists of ions of heavy gases and acts on the surface of the implant by removing any contaminants.


In fact, during treatment with argon (an inert gas that does not react with the titanium surface), the partially ionised atoms of the gas bombard and "sweep" the surface of the implant violently, removing the organic contaminants and thus obtaining a surface free from further traces and residues.

Another important step is represented by the class "C" cleanroom, according to GMP, which ensures the packaging of the dental implants in an environment where the concentration and quality of the particles in the air are kept constantly at the levels pre-set by the regulations in force in the medical field.

STERILISATION

The sterilisation process is the final step of the decontamination treatment and is carried out by beta ray irradiation. Sterilisation has a Quality Management System in compliance with UNI EN ISO 13485 standards and provides a guarantee of the product's sterility (certified duration of 5 years).

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The conical coupling guarantees a perfect antibacterial seal and high resistance of the connection between the components.

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SEM 60x images of the adhesion between the implant and the abutment.

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FEM analysis has proven the high resistance of the implant connection and the reduced transfer of non-homogeneous axial and lateral forces from the abutment to the implant and in turn to the bone, thereby increasing the likelihood of a successful treatment by reducing micro-movements and gaps between bone and implant.

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① *The conical implant-abutment coupling guarantees a better protection of the interface from bacterial and mycotic colonisation, by reducing the proliferation of bacteria.*

It also restricts or prevents the micro-movement of the abutment in relation to the implant when loaded, with consequent mechanical and biological advantages.

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② *Platform switching allows a better transfer of the load between implant, abutment and bone and the subcrestal positioning allows the bone to grow over the platform and therefore results in an implant system with a better biomechanical performance.*

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③ *The concave profile favours a greater quantity of soft tissue, better protection of the system and a lower risk of retraction, by constituting an effectively integrated biological system.*

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BIOLOY® Surface: Titanium surface micro structured by grit-blasting with a calcium magnesium carbonate $[\text{CaMg}(\text{CO}_3)_2]$ based media and then coated by a thin TiNbN layer.

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
SEM image at higher magnification of the microtopography of the BIOLOY® surface: macro-particles, or droplets, fused together are visible on the coating surface.

The coating modifies the typical sandblasted surface topography by rounding picks and valleys and by the presence of droplets.

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PVD is a process by which positively charged metal ions are evaporated in a vacuum chamber and react with nitrogen gas introduced to the chamber. The surfaces that are to be coated are negatively charged to allow a strong bond to form between the substrate (implants) and the coating.

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BIOLOY® Surface: cells seeded onto BIOLOY® Surface took a morphology characterized by both shapes seen on the machined titanium and on the sandblasted titanium surfaces. Cells spread over the smoother peaks and droplets of the BIOLOY® Surface.

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*p-value<0.05; **p-value<0.001.

Cell viability has been tested by AlamarBlue® assay.

Cell viability of cells cultured on BIOLOY® Surface is statistically higher than cell viability of cells cultured on both sandblasted titanium and machined titanium at 48 and 72 hours:

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Three histologic evaluation of samples from patients treated with a custom manufactured titanium mini-implants with BIOLOY® Surface inserted in the mandibular arch and removed at 3-4 months.

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BIC (bone to implant to contact): $46,8\% \pm 9,15$

BI (bone ingrowth): $59,8\% \pm 4,23$

Osseointegration of Titanium dental implants with BIOLOY® Surface:

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Image of the sample stained with Toluidine Blue and Pyronin G, which shows calcified bone in different phases of mineralization with different shades of blue/purple. Regenerated highly vital bone during its physiological remodelling can be observed within the coils. Total magnification 100x.

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Image taken at polarized light of the sample stained with Toluidine Blue which shows calcified bone in different phases of mineralization. Regenerated highly vital bone during its physiological remodelling can be observed within the coils as can be noted by the birefringency (polarized light), total magnification 100x.

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In vitro and in vivo study were performed at Department of Biomedical, Surgical and Dental Sciences, Università degli Studi di Milano, Milan, Italy.

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
CONNECTION IMPLANT - PROSTHETIC PLATFORM

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ROSSO=small

VERDE=regular e regular plus.

2 connection platforms/3 emerging profiles:

SMALL for Ø3,3 implants


REGULAR and REGULAR PLUS for Ø3,85/4,3/5/6 implants

to be chosen according to the space available and to the characteristics and quantity of the soft tissues.

A correct coupling between the implant and prosthetic part is ensured by a univocal colour code that identifies the appropriate platform for each implant. Each pack is colour-coded to match the platform:

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The purpose of these components is to prepare the gingival site to fit perfectly around the permanent prosthetic device. The prosthetic parts have the same profile as their corresponding healing screws.

		⊙ SMALL	⊙ REGULAR	⊙ REGULAR PLUS	
Vite di guarigione Conica <i>Conical Transmucosal Screw</i>	H.1,5	EKH11	EKH21	EKH31	
	 <i>Small Regular Regular Plus</i>	H.3	EKH13	EKH23	EKH33
		H.5	-	EKH25	EKH35
Vite di guarigione Cilindrica <i>Cylindrical Transmucosal Screw</i>	H.1,5	ECH11	ECH21	-	
	 <i>Small Regular Regular Plus</i>	H.3	ECH13	ECH23	-
		H.5	-	ECH25	-
Vite di copertura <i>Cover Screw</i>		⊙ SMALL ECC01	⊙ REGULAR EKC01		



Si consiglia un serraggio definitivo a 20 N-cm.
A final tightening at 20 N-cm is recommended.

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* During this phase, through different techniques, it is suggested that you should apply soft material that simulates the gum to facilitate the removal of the prosthesis part of the hardened mould, which will allow you to reproduce the soft tissue profile.

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**Peek Abutment has aesthetic function and must not remain in the oral cavity for more than 180 days.*

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These pre-formed components guarantee absolute precision and quality of mechanical and biological response. For better results, during the trimming phase, do not remove the titanium nitride coating from the transgingival part, because it is very useful and helps maintain the aesthetic result of the subgingival prosthetic protrusion.

The prosthetic parts have the same profile as their corresponding healing screws.

NOTE

An ideal gingival site is obtained using conical healing screws.

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The personalised abutments offer a number of clinical advantages. The bases for crowns are ideal for use with CAD/CAM technology, putting the dental technician in a position to plan and make the restoration with absolute precision, according to the dentist's instructions. The particular shape of the protruding profile guarantees precision when calculating the space needed to cement the crown. The trimmable abutments constitute an alternative to pre-shaped straight and angled abutments and leave the laboratory complete autonomy when trimming the device.

The prosthetic parts have the same profile as their corresponding healing screws.

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NOTE

The ideal gingival site is obtained through the use of conical healing screws.

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A final fixing of the prosthesis on the abutment of 15 N·cm max, straight CPA 30 N·cm and angled CPA 20 N·cm max.

REGULAR		H.1,5	H.3	H.4	H.5
Moncone CPA Dritto <i>CPA Straight Abutment</i>		E8421	E8423	-	E8425
Moncone CPA Angolato 17° <i>CPA Angled Abutment 17°</i>		E841721	E841723	-	-
Moncone CPA Angolato 30° <i>CPA Angled Abutment 30°</i>		-	E843023	E843024	-
Analogo <i>Analogue</i>	Transfer	Cuffia di protezione <i>Protective Guard</i>	Cilindro provvisorio in titanio <i>Temporary Cylinder in Titanium</i>	Cilindro calcinabile <i>Calcifiable Cylinder</i>	Vite M1,4 <i>Screw M1,4</i>
E8A2	E8T2	E8VP2	E8CTT2	E8CC2	E8V

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NOTE
The ideal gingival site is obtained through the use of cylindrical healing screws.

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1

This methodology entails resolving the total lack of teeth or other situations in which the teeth are compromised, inserting only 4 implants (in some cases there may be 6 implants) in the anterior region of the jawbones and/or mandible. The All on Four technique allows for the placement of a prosthesis for immediate loading called Toronto with a minimum of 10 teeth and then, after six months, to rehabilitate the patient with a definitive prosthesis, with a titanium structure, of 12 teeth.

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
Insert the CPA abutments into the implant by using the dynamometric ratchet and the tightening key, ensuring that the prosthetic axis is parallel. Tighten the fixing screw of the CPA Abutment manually with the screwdriver. Carry out the final tightening with the torque wrench.

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NOTE


Once all the straight and angled CPAs have been inserted, it is recommended that an X-ray should be performed to verify that the implants are properly connected.

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Screw a Transfer for CPA on each abutment and take the impression with a single open tray.

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Drill holes in the impression-tray according to the individual restoration so that the positioning screw of the impression-tray is protruding. Take the impression.

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
The dental laboratory can create the mould by using the CPA-Analogues that perfectly reproduce the conical head of CPA abutments.

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
Unscrew the Transfers from the pillar. During the temporary laboratory phases, it is possible to place the protective guard to maintain the soft tissues.

After the removal of the guard, it is necessary to verify that the CPA abutments are properly connected with the implants by the tightening with the Dynamometric Ratchet.

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
The temporary prosthesis can be made by using the temporary Cylinders in CPA titanium.

The CPA cylinders must not be cut below the first ring from the bottom. Tightening with the Dynamometric Ratchet at 20 N·cm.

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To make the permanent prosthesis, use the CPA Calcinable Cylinders. The permanent prosthesis must not be tightened with the torque higher than 20 Ncm.

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A final tightening of the prosthetic screw at 30-35 N·cm is recommended.

		⊙ SMALL	⊙ REGULAR
Monconi a Sfera <i>Spherical Abutment</i>	H.1,5	E4S11	E4S21
	H.3	E4S13	E4S23
	H.5	E4S15	E4S25

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The EQUATOR® system makes it possible to attach the prosthesis conveniently and securely to the patient through the use of elastic retention clips.

It is equipped with different retentive caps to be able to exert the desired retention force on the matrices.

Tightening the pillars must be carried out manually, with force, all the way up to stop, by repeating this operation two/three times to get a good adaptation of the male thread with the female component.

For the application of the prosthesis in the patient's mouth, please refer to the instructions for use of the manufacturer.

A final tightening of the prosthetic screw at 30-35 N-cm is recommended.

		⊙ SMALL	⊙ REGULAR
	Pilastro Equator® Equator® Pillar		
	H.1.5	ESQT01	ERQT01
	H.3	ESQT02	ERQT02
	H.5	ESQT03	ERQT03
	Cappette Ritentive Retentive Cap Kit		
	ALTA HIGH	VIOLA VIOLET	
	STANDARD	BIANCO WHITE	
	LEGGERA LIGHT	ROSA PINK	
	EXTRA-LEGGERA EXTRA-LIGHT	GIALLO YELLOW	
		Guscio in Acciaio Inox Stainless Steel Shell	

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NOTE

The ideal gingival site is obtained through the use of cylindrical healing screws.

KIT CHIRURGICO


SURGICAL KIT

Kit Chirurgico
Surgical Kit

KT-108

STRUMENTARIO

INSTRUMENTS

FRESE CHIRURGICHE

SURGICAL DRILLS

Fresa Lanceolata / Lanceolate Drill


	FIL01
---	-------

Numero di giri consigliato 800-1000 g/m-Fisiologica 140 ml/min.
 Numbers of turns recommended 800 -1000 g/m-Saline 140 ml/min.

Frese Cilindriche Cylindrical Drills

	CON STOP WITH STOP	SENZA STOP WITHOUT STOP
	FCE20	FCE20S
		
	FCE33	FCE33S
		
	FCE38	-
	FCE43	-
	FCE50	-
	FCE60	-

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Numbers of turns recommended 600 -800 g/m-Saline 140 ml/min.
 The stop-less drill is to be used in cases where there is little space, by using the laser notches as a reference for depth.

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It is used for the manually screwing the implants in combination with the extensions.

In case of uses where surgical residues are generated (blood, secretions, tissue residues), the key must always be removed, and then positioned in an adequate disinfectant solution.

This operation simplifies cleaning because dried residues are the cause of corrosion.

After cleaning, thoroughly rinse the parts in water and, by using a nylon brush, brush the inner and outer surfaces of the different parts of the key.

During the cleaning step, avoid any contact between the various parts of the key.

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Completely unscrew the adjustment ring of the control lever.

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Unscrew the screw of the fixing cover of the ratchet unit by using the adjustment ring that acts as a screwdriver.

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NOTE:
This tool is not a measuring device.

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RECOMMENDATIONS FOR THE TREATMENT OF DEVICES DURING AND AFTER SURGERY

Place the instruments used during surgery on a tray containing distilled water.

After surgery, the reusable devices must be rinsed to promptly remove the organic residues.

Wash all parts that make up the devices; disassemble the devices only where it is required.

The reusable devices must be thoroughly cleaned with a plastic toothbrush with stiff bristles that must not be made of metal. Use brushes and cleaners that have in turn been decontaminated, washed and sterilised.

Perform decontamination by immersing each device in a disinfectant solution suitable for the type of material. Rinse all devices with running water, possibly demineralised, to eliminate residual traces of detergent.

Dry them, which is necessary not to compromise the sterilisation process, with a clean soft cloth or with compressed and filtered air.

Immediately prior to sterilisation, the devices must be packaged with a suitable material.

Saturated steam autoclaving sterilising:

The devices, thoroughly dried after cleaning, must be suitably packaged and autoclaved according to the validated sterilisation process and according to the instructions provided by the manufacturer of the autoclave.

The effectiveness of saturated-steam autoclaving sterilisation was verified at 121°C at 1 pressure atmosphere.

Preparation of sterile devices:

Pack the individual instruments.

Put the instruments in the appropriate tray.

Pack the surgical tray.

Add a label with an indication of sterility.

Add a label with the expiration of sterility.

Recommendations

Fully immerse the surgical trays in the decontaminating solutions and subsequently in the detergent solutions.

Use only detergents at PH 5-9 (higher or lower PH values destroy the layer of anodised aluminium and any writings or symbols).

Always verify the correct functioning of the equipment for sterilisation.

Storage

After sterilisation, the devices must remain in the wraps used for sterilisation and can only be opened immediately before use.

The storage period of the sterilised items inside the wraps must not exceed the expiration of sterility.

WARNINGS

Never let the organic residues dry.

Do not put instruments of different metals in the ultrasound device.

Do not let sharp tools come into contact.

Do not use metal brushes.

After disinfection, rinse the instruments very thoroughly with running water.

Inspect the instruments that are not damaged and test their functionality.

Devices that are not perfectly clean may contain residues left on the surface of the instruments and make sterilisation ineffective.

The parts provided in the NON-STERILE package, including also those that transit for a short amount of time in the oral cavity, must be disinfected by using specific products for medical devices and STERILISED (autoclaved) consistently with the instructions for use provided by the manufacturer of the device.

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The background colour of the label identifies the colour code assigned to each coupling platform.

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The toning label (adhesive sticker in a material susceptible to the nuclear rays of the sterilisation cycle) is placed to ensure that irradiation is verified under all conditions.

Parte Protesica Prosthetic Part

REF. E421
EXACTIVE
REGULAR H.1.5
Ti6Al4V ISO5832/3
RIVESTIMENTO TIN

MONCONE DRITTO IN TITANIO
STRAIGHT ABUTMENT IN TITANIUM

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IMPORTANT: to avoid accidental dropping or loss of the implant, hold the container vertical at all times with the titanium part underneath and the transparent vial on top. The same position must be used during opening and implant removal.

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
**Lock the implant by pressing the support to correctly orient the internal shape of the implant with the outer shape of the inserter.

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PROCEDURE FOR INSERTING THE IMPLANT

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
Open the primary package and take out the sterile wrap. Once you have opened it, drop the container of the implant onto the sterile sheet of the operating table.

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Remove the implant directly from the package with the contra-angle screwdriver. The contra-angle inserter can be used as a manual inserter when the DGM article is used with the AD adaptor.

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Before inserting the implant, wash the receiving site and let it bleed. Inserting the implant can be done by using the contra-angle screwdriver.

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
Insert in the predetermined position. It is possible to use the dynamometric manual ratchet to finalise the position of the implant. If a higher torque is required, after using the inserter, use the contra-angle or ratchet driver up to a maximum torque of 45N-cm. Beyond these values, the tap should be used to facilitate insertion whilst observing the loading parameters for which the implant was designed.

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Before its positioning, wash the head of the implant and the surrounding areas with saline solution, and aspirate thoroughly. Then position the cover screw and tighten it manually with the appropriate driver.

Alternatively, place the healing screw or any other prosthetic restoration chosen by using the specific clamping force.

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ATTENTION: If you exceed the screwing threshold of 45 N-cm, excessive pressure may cause apical ischemia and lack of nourishment. In addition, the mechanical hexagonal male-female components could become deformed irreversibly.

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DRILLS

The drills are manufactured using stainless steel and require copious external irrigation during use.

Use with a reciprocating motor with input and output from 1 to 2 seconds.

Extract the drill without stopping the motor, in order to rinse any drilling detritus out of the implant site.

Irrigate the site to remove any bone fragments before proceeding with the further drilling.

Drill as far as the desired depth, which will depend on bone quality and type and the diameter of the implant.

Drills may be up to 1 mm longer than the implant and corresponding depth notches.

This length should be considered when drilling close to vital anatomical structures.

COUNTERSINK

The countersink can be used to adapt the housing for the head of the implant to the receiving site, (bearing in mind bone quality and thickness at all times).

In the presence of soft bone or when maximum primary stability is required, use of the Countersink can be avoided by respecting the torque values indicated.

THREADING TAP


In the presence of medium or dense bone, the tap must be used to create a thread, in order to enable the insertion of the implant in the site prepared.

When using the tap, apply the necessary pressure during insertion to reach the depth required, then invert the direction of rotation for extraction.

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permedica s.p.a. has a Quality Management System that complies with the requirements of the European Directive 93/42/EEC (as subsequently amended and supplemented) relating to medical devices and certified by the certifying agency Italcert according to international standard ISO 13485.

The Quality Management System pays particular attention to the application of the directive and the specific regulations on medical devices, to guarantee absolute control over all phases of planning, production and sale.

Our processes are subject to continuous updates and innovations, for the purpose of constant improvement. This catalogue provides users with a guide to optimise the preparation of an implant site and the insertion of the implants of the EXACTIVE® line.

This must not to be understood as an alternative to the doctor's training and professional experience. This catalogue incorporates, without replacing it, the Instructions for Use sheet that accompanies each permedica s.p.a. device.

Each device is identified and can be ordered by using the item code listed in this catalogue.

For further information or clarifications please contact your local dealer or manufacturer.

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The label of each device contains data relevant to the types of materials/coatings used.

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The materials used for implantable devices produced by permedica are the following:

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The materials used for the prosthesis part are the following:

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- > *Alloy Ti 6Al4V (ISO standard 5832/3)*
- > *Alloy Ti 6Al4V (ISO standard 5832/3)*
- + *TiN coating*
- > *PEEK (polyetheretherketone)*
- > *PMMA (polymethylmethacrylate)*

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The success of permedica's implantology systems is the result of extensive studies on design and production technologies.

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deutschland GmbH
permedica deutschland GmbH