zubler CONCEPT VEST™ LiSi_2

Speed Investment material for Lithium Disilicate Press ceramics

MIXING RATIO FOR CROWNS, VENEERS, PARTIAL CROWNS AND PRESS OVER TECHNIQUE: 80% LIQUID


	powder	liquid	dest. water	reaction time
Small	100g	20,8ml	5,2ml	30 minutes
Large	200g	41,6ml	10,4ml	30 minutes

MIXING RATIO FOR INLAYS: 50% LIQUID (Siliconring)

	powder	liquid	dest. water	reaction time
Small	100g	13,0ml	13,0ml	30 minutes
Large	200g	26,0ml	26,0ml	30 minutes

1. Fill liquid/water-mixture and investment powder both into a mixing bowl.
2. Start timer when bringing together the powder and the liquid.
3. Mix the investment by hand with a spatula.
4. Put the investment under vacuum for 15 sec. without starting the mixer.
5. Mix for 60 sec. under vacuum by 420 r.p.m.. If the RPM is not adjustable with your mixing device, please increase the mixing time up to 90 seconds.
6. After the mixing procedure, leave the investment under vacuum for 5 sec. .
7. Fill the ring former during vibrations.
8. After filling the ring former, protect it from any movement.
9. Remove the ring former five minutes before end of reaction time.
10. Edges of the ring should be scraped and base should be levelled.
11. Put the ring with the opening downside into the preheated furnace at 850°C / 1562°F .
IMPORTANT: The furnace must be clean! Please place the ring as centric as possible into your burnout furnace or rather by placement of divers rings please keep a distance of 2,5 cm among the rings and to the heat insulation (heating elements) of your burnout oven.
12. If there is more than one ring, keep an interval of 10 min. before the next ring can be inserted.
13. The recommended preheating time for a 100g ring is minimum 45 min., for a 200g ring 60 min., measured from reaching final temperature (850°C); each additional ring increases the time by about 10 min.
14. At the end of the preheating time, start the procedure of pressing.

Process and store investment only at a room temperature of 20°C / 68°F.
Differences in temperature can cause changed working and bench set times and an increase of expansion.

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